## Patient's and Physician's Manual

U.S. federal law restricts this device to sale by or on the order of a physician.

Read  
instructions  
before use

For  
single  
use only

**ATOS**  
**MEDICAL**

## **Caution**

United States Federal law restricts devices identified herein with the symbol “Rx” to sale, distribution and use by or on the order of a licensed practitioner trained in post-laryngectomy vocal and pulmonary rehabilitation. The availability of this product without prescription outside the United States may vary from country to country.

## **Disclaimer**

Atos Medical makes no warranty, either expressed or implied, as to the lifetime of the product delivered to the purchaser hereunder, which may vary with individual use and biological conditions. Further Atos Medical makes no warranty of merchantability or fitness of the product for any particular purpose.

## **Patents and Copyrights**



Remove and Skin-Prep™ are registered trademarks of Smith & Nephew United Inc.,

Provox® is a registered trademark owned by Atos Medical AB, Sweden. Provox® Micron HME™, Provox® LaryButton™, Provox® LaryTube™, Provox® LaryClip™ and Provox® TubeHolder™ are trademarks owned by Atos Medical AB.

Provox® HME has US patents (5.738.095; 5.042.468), other patents and patents pending.

# Contents

ENGLISH	5
DEUTSCH	15
NEDERLANDS	25
FRANÇAIS	35
ITALIANO	45
ESPAÑOL	55
PORTUGUÊS	65
SVENSKA	75
DANSK	85
NORSK	95
SUOMI	105
ΕΛΛΗΝΙΚΑ	115
ORDERING INFORMATION	125

## Descriptive Information

### Indications for Use

The Provox® Micron HME™ is a specialized stoma cover that acts as a heat and moisture exchanger (HME) and air filtration device for laryngectomized patients. Provox Micron HME **partially restores lost breathing resistance due to laryngectomy.** For patients with a voice prosthesis or surgical speech fistula it may also facilitate voicing.

Provox Micron HME is intended to be used with the attachment devices in the Provox HME System (*fig 1*).

### Contraindications

This device is contraindicated for patients with decreased level of consciousness, or reduced mobility of the arms and/or hands, who are unable to remove the device.

*For more information about attachment devices, please refer to the section “Attachment devices”.*

*For ordering information please refer to the section “Order information” on the last page.*

### Description of the device

#### General

Provox Micron HME is a disposable device that features a calcium chloride-treated foam core in a plastic housing. Provox Micron HME also contains electrostatic filter material which is efficient in blocking airborne particles (see **technical data** below). The housing of Provox Micron HME can be pressed down to occlude the stoma in order to speak with a voice prosthesis. When the pressure is released, the housing automatically retracts and opens the airway passage.

## **Electrostatic filtration**

Provox Micron HME is not only an HME, Provox Micron HME is also equipped with a layer of electrostatic filtration material. The filter consists of synthetic fibres which, during the manufacturing process, generate electrical charges on the fibre surface. The result is a highly stable, coarse fibred electrostatic filter media with high efficiency against sub-micron particles (see Technical data) and low air flow resistance.

## **Technical Data**

### **Measurements:**

Height	20 mm
Length	43 mm
Width	53 mm

<b>Weight:</b>	4 g
----------------	-----

Pressure drop after 1 h at: (According to ISO 9360)	30 l/min 60 l/min 90 l/min	0.7 hPa 1.9 hPa 3.5 hPa
Pressure drop after 24 h at: (According to ISO 9360)	30 l/min 60 l/min 90 l/min	0.7 hPa 1.8 hPa 3.5 hPa

Moisture output: (According to ISO 9360)	22 mg/L
---	---------

Moisture loss at VT=500 ml: (According to ISO 9360)	22 mg/L
--	---------

Bacterial Filtration Efficiency (BFE): (According to ASTM F2101)	≥ 99.8%
---	---------

Virus Filtration Efficiency (VFE): (According to MIL-M-36954C and ASTM F2101)	≥ 99.8%
--	---------

## **Expectations of the device**

### **Air Filtration**

Provox Micron HME provides capacity to filter inhaled air thereby reducing inhaled particles through consistent normal use.

The Integral filter helps to filter small airborne particles, e.g. bacteria, viruses, dust and pollen from passing through the device.

### **HME Function**

Continuous use may help improve pulmonary function and respiratory problems, e.g., coughing and mucus production. New users may experience slight initial discomfort related to increased breathing resistance. During the first weeks of use, mucus production may appear increased due to thinning of the mucus with retained water.

### **WARNING**

Unintentional or accidental pressure on the housing of the Provox Micron HME may cause difficulty in breathing or suffocation.

Inform the patient, caretakers and others of this feature to ensure that they understand the closing function. By putting pressure on the housing of Provox Micron HME, the airway passage is closed. Closing the airway in order to allow voicing is a well known feature for the laryngectomized patient with a voice prosthesis, but may be unknown for patients without voice prosthesis.

## **PRECAUTIONS**

Each disposable Provox Micron HME device can be used for a maximum of 24 hours (after opening the package) in order to maintain humidification and filtration capacity. Using the product exceeding the specified 24 h intervals may cause colonization of harmful microorganisms in the product.

Provox Micron HME provides filtration of the breathing air through consistent normal use. A fresh device on an airtight seal gives good protection (*read section "Attachment devices" for instructions how to ensure an airtight seal*). However, since there are other pathways for pathogens to enter the human body Provox Micron HME **can never guarantee total protection**.

For effective humidification and filtration the breathing air must pass through the Provox Micron HME. Check that the seal is airtight by closing the Provox Micron HME (*fig 2*). Restore the seal(s) if there is indication of leakage.

**Do not reuse, clean or soak** Provox Micron HME with water or other liquids. The function of the heat and moisture exchanger will be substantially reduced if the calcium chloride is washed out of the foam, and the filter efficacy will be compromised.

**Do not disassemble** Provox Micron HME. Disassembly will destroy its function.

## **Instructions for use**

### **Insertion**

Provox Micron HME can easily be inserted and removed manually from the attachment devices in Provox HME System when needed (*fig 3*).

## **Using Provox Micron HME to close the tracheostoma for voicing**

In order to obtain an airtight seal for speech with a voice prosthesis or surgical speech fistula, press the housing down with a finger. This will direct the exhaled air through the voice prosthesis. Once the finger is lifted and the pressure is released the airflow will pass through Provox Micron HME again (*fig 2*).

You can also use the device to close the stoma if a pressure increase of the air in the lungs is needed (e.g. for coughing).

### **Coughing**

Provox Micron HME can be removed during or after coughing if the stoma needs to be cleaned from mucus. Generally, the attachment devices in Provox HME System can stay in place during cleaning of the stoma or during Cassette replacement.

### **Replacement**

Replace Provox Micron HME when needed (*fig 3*). To ensure its proper function, the same Cassette should not be used for **more than 24 hours** (after opening the package). Cassette replacement frequency is based on individual patient needs.

### **Storage**

Store the device in a clean, dry and dark area at room temperature (15°C–25°C/60F–80F).

### **Re-order information and instructions on accessories**

For order codes, please refer to the table on the last page.

## Other information

For information about clinical studies, on how to contact your local representatives, and to get information on other Atos Medical products, contact our customer service department or visit us at [www.atosmedical.com](http://www.atosmedical.com).


## Date of Printing

See version number on the back cover of this Manual. See version number on the back cover of this Manual.

## Attachment devices

The Provox Micron HME and other Provox HME Cassettes are intended to be used with the following attachment devices and accessories in Provox HME System.

For more information read Instructions for use for each of the products.

## **Provox® Adhesive**

Provox Adhesive base plates are single use devices that are intended to hold Provox Micron HME and other Provox HME Cassettes in front of the tracheostoma and to provide an airtight seal. The self-adhesive tape discs are supplied on a peel-away paper liner and come in two different shapes, i.e. oval (anatomical) and round, and three different types of glue: OptiDerm, FlexiDerm and Regular.

A special version, Provox XtraBase, is a Flexi-Derm type adhesive that has a conically shaped base. It is designed to facilitate and improve skin attachment for deep stomas and / or if an automatic tracheostoma valve Provox FreeHands HME is used.

OptiDerm (*fig 1a*) is intended for patients with a sensitive skin or if the skin is not yet healed properly, e.g. after surgery and/or irradiation. It absorbs moisture, is gentle to the skin and needs less frequent replacement. It can be applied gently on the 1st postoperative day. Average durability 36 hours.

FlexiDerm (*fig 1b*) is very flexible and has the strongest adhesive properties. It is especially suitable for patients with a deep and/or irregular stoma or who prefer a very soft and flexible adhesive. Due to the type of glue used, a suitable skin cleanser is recommended to clean the skin of any residual glue. Average durability 48 hours.

Regular (*fig 1c*) is a transparent, perforated adhesive. It is less flexible and sticky than FlexiDerm, but leaves less residual glue and is therefore easier and more comfortable to use. This adhesive is mainly intended for persons with a regular stoma and skin. Average durability 24 hours.

## **WARNINGS**

- The adhesive glue may irritate the skin. Stop using the adhesive if skin irritation develops and consult your physician.
- The Provox HME Adhesives should not be used during radiotherapy of the skin around the stoma or for 2 weeks after radiotherapy.
- Do not use the OptiDerm adhesive as attachment for the Provox ShowerAid. The glue may dissolve in contact with water and the Adhesive may become loose, allowing water to enter the trachea (*fig 4*).

## **Operating Information**

### **Set-up Instructions:**

1. Clean the stoma and the skin around it with Provox Cleaning Towel and/ or soap and water, and then dry the area carefully.
2. You may apply skin preparation agents, like Skin-Prep™, or, if needed, Provox Silicone Glue, according to their instruction for use.
3. If OptiDerm is used: Please carefully pre-warm the Adhesive base plate between your hands for improved attachment.
4. Remove the protective paper from the adhesive. Stretch the skin gently away from the stoma, and then apply the adhesive, in order to get optimal attachment to the skin, even in skin folds.
5. Massage the adhesive carefully gently for approximately one minute in order to improve skin attachment.

## **Replacement of the Adhesive base plate**

The Adhesive base plate can stay in place as long as it is glued properly to the skin. To facilitate removal of partially attached Adhesives, the use of adhesive remover, e.g. Remove™ may be helpful. After Remove™ is used the skin should always be cleaned with a cleaning towel or soap and water. Read instructions. Before applying a new Adhesive base plate, the skin should be carefully cleaned and dried.

**CAUTION:** When cleaning the skin from residual glue, prevent particles/fluids from entering the stoma. In case of using OptiDerm in the postoperative period or sensitive skin, the Adhesive should only be replaced when completely loose.

## **Storage**

Store in a clean, dry and dark area at room temperature (15°C–25°C/60F–80F).

## **Provox® LaryButton™**




Soft silicone laryngectomy stoma buttons (*fig 1d*).

- with or without using Provox LaryClip or Provox TubeHolder.  
Read separate instructions.

## **Provox® LaryTube™**

Soft silicone laryngectomy tubes (*fig 1e*).

- with or without using Provox LaryClip or Provox TubeHolder.  
Read separate instructions.

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© Atos Medical AB, 200807A, REF 10029.