FreeHands HME[®]









Instructions for use









Read
instructions
before use

ATOS

MEDICAL

ENGLISH..... 4-15

DEUTSCH 16-29

NEDERLANDS..... 30-41

FRANÇAIS..... 42-54

ITALIANO..... 55-67

ESPAÑOL 68-80

PORTUGUES..... 81-92

SVENSKA 93-103

DANSK..... 104-114

NORSK..... 115-125

SUOMI..... 126-136

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TÜRKÇE..... 152-162

Contents index

General Information	4-8
Clinicians.....	9-11
Patients.....	12-15

Provox® FreeHands HME® - Instructions for Use

General Instructions

These Instructions for Use contain information for the clinician and the user. Read this information carefully each time a Provox® FreeHands HME® is provided. There may be new information. This information is not intended to take the place of discussions between clinician and user. Make sure any questions you may have about the use of the Provox®, FreeHands HME® are answered.

Before Use

Before the device can be used **for the first time**, the proper function of the cough-relief valve must be checked and adjusted if necessary, and the right speech valve membrane chosen **by a clinician** trained in voice and pulmonary rehabilitation. Readjustment to individual conditions may be needed from time to time. A check-up after some weeks of use is recommendable.

Prior to each use, inspect the product to ensure it is undamaged and that all parts and accessories are complete and in proper condition. Do not use the product if any components show signs of cracking or other structural damage.

Product description

The product consists of two functional units that work together:

- a Provox® FreeHands speech valve unit in a reusable housing, containing two independent valves: an internal, exchangeable valve membrane (speaking valve) and a hinged valve on the top (cough-relief valve)
- a disposable Provox® FreeHands HME® Cassette, which has to be connected to the lower part of the valve housing before it is ready to use.

After proper adjustment and attachment of speech valve unit and HME cassette, the completed Provox® FreeHands HME® is inserted into either Provox® Adhesive base plate or Provox® LaryTube with ring. (The latter is not included in the package.) Accessories included with the Provox® FreeHands HME® are:

- five disposable Provox® FreeHands HME® Cassettes
- two Provox® XtraBase® Adhesives base plates
- a safety screwdriver for the adjustment of the cough-relief valve
- a plastic forceps for replacement of the speech valve membrane
- two extra speech valve membranes: "light" (white dot) and "strong" (green dot). The "medium" membrane (blue dot) is already mounted in the valve.
- a plastic cleaning and storage box for the reusable valve unit
- one package of Skin-Prep™ wipes
- Instructions for Use

How does the Provox® FreeHands HME® work?

The speech valve unit allows hands-free tracheostoma occlusion for voice prosthesis users.

The speech valve closes as exhalation pressure is increased briefly, just as it increases when one starts to talk loudly in normal laryngeal speech. If the exhalation pressure does not increase rapidly, the valve stays open and allows normal breathing.

Once the valve is closed, a small magnet holds it in the closed position. There is no need to maintain high air pressure in the trachea/windpipe in order to keep the valve closed.

When the pressure in the trachea/windpipe drops to a very low level, or during inhalation, the valve opens and allows free breathing again.

There are two valve positions, which can be adjusted by simple rotation of the valve housing. The "speaking position" (On) allows the automatic change between stoma closure and normal breathing as described above. The "breathing position" (Off) allows a higher airflow without closing the valve. This is also achieved by the use of two magnets, which keep the valve open despite the increased airflow during physical activity.

The cough-relief valve helps reduce the risk of trauma to the airways, which might occur due to sudden pressure increases while coughing. The cough-relief valve also helps to diminish the pressure on the Adhesive and, thus, premature loosening of the airtight stoma seal. The pressure (cough intensity) needed to

open the cough-relief valve can be adjusted individually with the included screwdriver before the Provox® FreeHands HME® is put in place.

The HME unit of the Provox® FreeHands HME® is a full heat-and-moisture-exchanger, which is used for partial restoration of the lost upper airway functions. It also separates the speech valve from direct contact with the airway and, thus, prevents mucous from being coughed into the valve part of the device. It also reduces the risk of inhalation of small particles, e.g. coming from a defective valve.

Technical Data

Pressure drop at 30 l/min:

81 Pa (0,83 cmH₂O)

Pressure drop at 60 l/min:

242 Pa (2,47 cmH₂O)

Moisture Loss: at 500 ml tidal volume

19,1 mg/l

(according to ISO 9360-2; 2001)

Indications

The Provox® FreeHands HME® is an accessory to prosthetic vocal and pulmonary rehabilitation after a total laryngectomy. It is designed to allow automatic tracheostoma closure (instead of digital stoma occlusion) for voice prosthesis users (e.g. Provox®, Voice Prosthesis), while providing heat-and-moisture-exchange (HME) for pulmonary rehabilitation.

The device should be adjusted for the end user by a clinician trained in voice and pulmonary rehabilitation, e.g. a speech pathologist, before it can be used without supervision.

THE PROVOX® FREEHANDS HME® IS INTENDED FOR SINGLE PATIENT USE.

Contraindications

The device must not be used if there is no properly functioning voice prosthesis in place or if the laryngectomee is unable to speak with such a device. The user must be able to produce intelligible voicing by means of digital stoma occlusion. Consequently, all contraindications for prosthetic voice rehabilitation also contraindicate the use of a Provox® FreeHands HME®. Absolute and relative contraindications for prosthetic voice rehabilitation include, but are not limited to: Radiotherapy doses greater than 70 Gy in 7 weeks, poor overall physical condition, heart insufficiency, coagulopathy or anti-coagulant therapy, mental diseases, hypersensitivity to the material of the voice prosthesis.

GENERAL

Allergy or hypersensitivity against any of the following materials contraindicates the use of the Provox® FreeHands HME®: metals, polycarbonate, compound thermoplastics and silicone.

Lung emphysema, asthma and other pulmonary diseases, which affect lung capacity, are contraindications for breath-controlled automatic speech valves. Patients suffering from such diseases should only use the Provox® FreeHands HME® if it has been determined by their physician that the benefits of the device clearly outweigh the risk of its use.

The stoma size and the condition of the peristomal skin may be factors contraindicating the use of the device and its required accessories. A stoma, which is too small can limit the inhaled airflow and thereby lead to too low a tidal volume. Hypersensitive skin may make it impossible to attach the valve properly.

Warnings

Always assemble the device (HME cassette and valve unit) before use. Using not properly assembled or single device components may increase the risk of aspiration of these components or device fragments, e.g. of lose magnets in the case of a speech valve unit defect. Using the Provox® FreeHands HME® Cassette without the valve unit may lead to aspiration of the cassette or parts of it.

If the connection between the speech valve unit and HME cassette becomes loose during use or removal of the device, the HME cassette must be removed from the tracheostoma only by removing the complete Adhesive base plate from the skin. No other manipulation must be done in such situation, in order not to increase any risk of aspiration.

Excessive coughing-up of mucus may suddenly block or nearly block your airway. Make sure you always have the possibility to remove the device quickly from your tracheostoma if needed.

Do NOT remove the foam from the Provox® Free-Hands HME® Cassette. Removal of the foam eliminates the physiologically needed HME features, damages the valve due to direct mucus contact and increases risks for the user, including the risk of incidental loss of small device parts and the risk of aspiration.

Do NOT reuse the foam, e.g. by washing out. Doing so eliminates the calcium chloride needed for a sufficient humidity exchange and can increase the risk of infection.

Do NOT wear the Provox® FreeHands HME® while sleeping. Instead, use a Provox® HME Cassette with a digitally operated speech valve (picture P19) or use the

GENERAL

accessory Provox® HME Cap (REF 7730) with a Provox® FreeHands HME® Cassette (see Accessories).

The device is NOT MRI-compatible and must NOT be worn during radiotherapy or x-ray.

Precautions

Handle the device with care.

Use only the recommended cleaning solutions/tabs-lets for cleaning and care. We recommend commercially available denture cleaning tablets.

Never use oil-based detergents, boil or sterilize the device, since this might destroy the valve.

For adjustment of the valve, i.e. the exchanging the speech valve membrane or adjusting the cough valve magnet, **only use the Adjustment Kit** (plastic forceps and safety screwdriver) provided with the Provox® FreeHands HME®.

Adverse Events

The patient should watch for any change in the stoma tissue, the tissue around the voice prosthesis, increasing difficulties in breathing, or similar symptoms. If any of these symptoms or any other difficulties possibly related to the use of the device occur, a clinician should be consulted.

Some patients may experience a reaction to the Provox® Adhesive. Discontinue use of Adhesive until the reaction resolves. The use of Provox® OptiDerm (hydrocolloid) Adhesive, skin protection with Skin-Prep™ or wearing a LaryTube may reduce the risk of such reactions.

Device Lifetime Expectation

The HME should be replaced at least every 24 hours or as often as needed when the HME is blocked by mucous, etc.

Device ageing and long-term usage of the reusable valve unit may lead to improper valve function. Symptoms for a worn-out Speech Valve Membrane can include delayed closing and / or opening, air leaks and increased closing sound. Symptoms for a defective valve housing include mechanical defects and / or changes in the cough-relief valve function.

Regular conditions provided, the Speech Valve Membranes withstands a minimum of 6 months use, and the valve housing of 12 months use. Proper handling, cleaning and storage according to the instructions further below are important factors to maintain the device function.

Speech Valve Membranes for replacement are separately available (REF 7713, 7714 and 7715).

Atos Medical makes no warranty as to the lifetime of the device, which may vary with individual use and biological conditions.

Instructions for the Clinician

See also the illustrations for clinicians at the back of the manual **[C1 through C17]**.

If you are a patient, please refer to “Instructions for the Patient” further below.

Adjustment

STEP 1 – PREPARATION OF THE STOMA

Carefully clean and dry the skin around the stoma using Provox® Cleaning Towel or water and a mild soap. Make sure that the skin is dry before proceeding [picture C1]

Apply Skin-Prep™ to the skin around the stoma before you put on the Adhesive to enhance the sticking properties of the Provox® Adhesive. Let it dry for about 30 seconds and then apply the Adhesive base plate. [picture C2]

Provox® Silicone Glue (REF 7720), can be used to prepare the skin before applying the Adhesive base plate, thus providing a better seal. (The silicone glue is not included in the Provox® FreeHands HME® package.)

STEP 2 – PREPARATION OF THE DEVICE

Place a Provox® FreeHands HME® Cassette on a flat surface, with the plastic grid down and the open foam surface on top. Place the valve housing on top of the HME Cassette and push it down until it clicks, and the Cassette is tightly fixed to the housing. Check the proper fit manually. [picture C3]

Hold the device with the HME Cassette facing away from you. When looking at the top of the speech valve unit, twist the housing counter-clockwise until it stops. The valve is now in the “breathing position” (OFF). [picture C4]

STEP 3 – INSERTION AND ADHESIVE SEAL

Insert the Provox® FreeHands HME® into the Provox® Adhesive base plate and ask the patient to breathe normally. [picture C5]

Inspect the Adhesive seal. There should be no leakage. This can be achieved by massaging the Adhesive base plate onto the skin. If there is still leakage, a new Adhesive base plate should be used and carefully ap-

plied to the skin. Skin-Prep™ and silicone glue can be used to enhance the seal of the Adhesive base plate. **[picture C6]**

STEP 4 – SELECTION OF THE SPEECH VALVE MEMBRANE

Turn the valve “ON” to the “speaking position” (clockwise, if the patient is sitting opposite you), **[picture C7]**.

After a short time of adaptation and maybe some attempts, starting with low pressure, the patient should be able to speak “hands-free” without too much effort, **[picture C8]**. If this is achieved, proceed to step 6.

The valve originally comes with the medium (blue dot) membrane in place, to allow in many cases for immediate use without further adaptation.

If this fails, consider the following possibilities:

A) The valve closes too late, too briefly or not properly, allowing air to leak along the valve. The voice is weak and speaking is possible only for a very short time and much effort: Go to step 5 and chose a valve membrane with higher flexibility, which is easier to close (light or medium, respectively).

B) The valve closes too easily, e.g., it sometimes closes as a result of deep breathing. Go to step 5 and chose a valve membrane with lower flexibility (strong or medium, respectively).

C) The cough-relief valve opens during speech as might happen with patients used to other, conventional valve systems. Go to step 6, adjust the cough-relief valve, and repeat step 4.

STEP 5 – REPLACEMENT OF THE SPEECH VALVE MEMBRANE

If the membrane has to be replaced, use only the plastic forceps provided with the device for manipulation, in order to avoid damaging the membrane. **[picture C9 and C10]**

Lift the cough-relief valve (lid) with your fingers and hold it open. **[picture C11]**

Grasp the membrane next to the fixation point and pull it out. **[picture C12]**

Hold the cough-relief valve open and insert the new membrane by grasping it next to the fixation point and pushing it in. **[picture C13]**

STEP 6 – ADJUSTING THE COUGH-RELIEF VALVE

Once the speech valve membrane has been chosen, the patient should be asked to cough with increasing force. The valve should open as soon as the pressure increases to a level of discomfort during the cough.

Properly adjusted, it should not open during loud speech. During both speaking and coughing, the cough-relief valve should allow the Provox® Adhesive to remain in an airtight position. [picture C14]

The default setting of the cough-relief valve is at a medium opening force. However, because of individual differences, it is adjustable in order to keep the valve closed even during loud speech and to permit the valve to open instantly during a sudden cough.

To adjust the cough-relief valve, remove the Provox® FreeHands HME® from the Adhesive base plate and hold it with the HME part facing to you. [picture C15]

Locate the position of the cough-relief valve magnet (below the point on the top of the cough-relief valve) and, using only the safety screwdriver provided with the device, push the HME foam to the side, enough to see and reach the magnet screw in the valve housing. [picture C16]

The magnetic force of the valve can now be made stronger (by turning the screw clockwise) or weaker (by turning the screw counter-clockwise). Start with an adjustment of about 90°. [picture C17]

Stop immediately if resistance is felt while turning the screw, otherwise you can damage the valve irreparably.

STEP 7 – INSTRUCTING THE PATIENT

The patient should be instructed on the proper use and daily care of this device as described below. This manual must be given to the patient.

Cleaning the Device

REMOVAL OF THE PROVOX® FREEHANDS HME CASSETTE

See step 5 in the Instructions for the Patient

NORMAL CLEANING OF THE REUSABLE SPEECH VALVE UNIT

See step 6 in the Instructions for the Patient

Disinfecting the Device

The Provox® FreeHands HME® is intended for single patient use. However, should disinfection become necessary, use the following procedure:

Cleaning prior to disinfection

Use the cleaning container, which comes with the Provox® FreeHandsHME®.

Remove the speech valve membrane from the valve unit.

Fit the valve unit onto the holder, the lid must be open.

Fill the cleaning container with a solution of luke-warm tap water and a commercially available denture cleanser tablet.

Immerse the membrane and the valve (attached to the holder) for 10 to 15 minutes (read instructions of the denture cleanser). The container may be shaken to facilitate cleaning.

Rinse the valve, membrane and container carefully under tap water.

Shake off the water.

Chemical Disinfection

Fill the container with fresh 70% Ethanol, minimum temperature 18° C / 65° F.

Immerse the valve (with open lid) and membrane. Shake the container for a short time and let all soak for a minimum of 30 minutes.

Tip out the Ethanol and let the valve, membrane and container air-dry.

Storage

Put valve and membrane into the cleaning container and leave that closed until the next use.

Instructions for the Patient

See also the illustrations for patients at the front of the manual [P1 through P18].

This patient section of the Instructions for Use describes the more regular (every-day) use of the Provox® FreeHands HME®. For all general information, e.g., Product Description, Indications, Contraindications, Warnings and Adverse Events, Device Lifetime Expectation, please review the general information section at the beginning of this manual.

Please make sure that you are fully instructed and that your clinician has adjusted the device, prior to starting to use the Provox® FreeHands HME®. In case of any problems possibly related to the use of the Provox® FreeHands HME®, please consult your clinician.

STEP 1 – PREPARATION OF THE STOMA

Carefully clean and dry the skin around the stoma using Provox® Cleaning Towel or water and a mild soap. Make sure that the skin is dry before proceeding [picture P1]

Apply Skin-Prep™ to the skin around the stoma

PATIENTS

before you put on the Adhesive to enhance the sticking properties of the Provox® Adhesive. Let it dry for about 30 seconds and then apply the Adhesive base plate. [picture P2]

Provox® Silicone Glue (REF 7720), can be used to prepare the skin before applying the Adhesive base plate, thus providing a better seal.

STEP 2 – PREPARATION OF THE DEVICE

Place a Provox® FreeHands HME® Cassette on a flat surface, with the plastic grid down and the open foam surface on top. Place the valve housing on top of the HME Cassette and push it down until it clicks, and the Cassette is tightly fixed to the housing. Check the proper fit manually. [picture P3]

Hold the device with the HME Cassette facing away from you. When looking at the top of the speech valve unit, twist the housing counter-clockwise until it stops. The valve is now in the “breathing position” (OFF). [picture P4]

STEP 3 – INSERTION AND ADHESIVE SEAL

Insert the Provox® FreeHands HME® into the Provox® Adhesive base plate (REF 7251-56, 7265) and breathe normally. [picture P5]

Inspect the Adhesive seal with means of a mirror. There should be no leakage. This can be achieved by massaging the Adhesive base plate onto the skin. [picture P6] If there is still leakage, a new Adhesive base plate should be used and carefully applied to the skin. Skin-Prep™ and Provox® Silicone Glue can be used to enhance the seal of the Adhesive base plate. (The silicone glue is not included in the package.)

STEP 4 – USING THE VALVE

ON (= automatic function; i.e., the valve lets you breathe and talk)

To switch the valve ON, twist the housing of the device in the opposite direction as in step 2, until you feel it stop, [picture P7]. You should now try to speak some words, beginning with low volume and pressure. In case you feel increased breathing resistance after coughing, you may have to take out the whole device for wiping the mucous away, or to replace the HME cassette with a new one.

OFF (= only breathing possible; for use e.g. during harder physical activities)

For heavier breathing turn the valve into the OFF position (like in step 2, you can do that while the valve is attached to the stoma). [picture P8]. This should avoid unwanted valve closure and also minimize the breathing

resistance as much as possible.

STEP 5 – REMOVAL OF THE PROVOX® FREEHANDS HME® CASSETTE

The HME cassette should be replaced at least every 24 hours or as often as needed when the HME is blocked by mucous, etc.

Take the Provox® FreeHands HME® out of the Adhesive base plate. [picture P9]

Hold the housing with one hand and compress the HME part with a little force, at a place where no crossbar touches the outer rim. [picture P10] The Provox® FreeHands HME® Cassette will be deformed during this step and come loose from the valve housing. [picture P11] The cassette must be thrown away. [picture P12] But the reusable speech valve unit is cleaned and prepared for its next use (step 6).

STEP 6 – CLEANING THE REUSABLE SPEECH VALVE UNIT

Use the cleaning container, which comes with the Provox® FreeHands HME®. [picture P13]

Fit the valve unit onto the holder, the lid must be open. [pictures P14 and P 15]

Fill the cleaning container with a solution of luke-warm tap water and a commercially available denture cleanser tablet.

Immerse the valve (attached to the holder) for 10 to 15 minutes (read instructions of the denture cleanser). [picture P16] The container may be shaken to facilitate cleaning.

Rinse the valve and container carefully under tap water. [picture P17]

Shake off the water.

Leave the valve in the holder of the container lid and let it air dry before the next use. [picture P18]

Caution: Never use oil-based detergents, do not boil or sterilize the device, since this might destroy the valve.

Replacement of the speech valve membrane

If the membrane has to be replaced, please refer to step 5 in “Instructions for Clinicians”. In case you’d like to change to another type of membrane, please consult your clinician to adjust the valve properly.

Adjustment of the cough-relief valve

In order to ensure a proper valve function, your clinician should do this adjustment.

Accessories: Provox® HME Cap (REF 7730)

The Provox HME Cap is a dome-shape titanium

PATIENTS

ring that allows you to use the special HME cassettes of Provox FreeHands HME® (REF 7712) without a FreeHands valve unit. This is especially of advantage, if you want to continue HME use during the night but have no other HME cassettes available than FreeHands HME® cassettes. Read instructions.

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7718 Provox® FreeHands Cleaning and Storage Box (1 pc)

7719 Provox® FreeHands HME Adjustment Kit (1 kit)

7720 Provox® Silicone Glue (1 container)

7722 Rx Provox® FreeHands HME Replacement Device
(Valve only)

7716 Rx Provox® FreeHands HME Replacement Device
(Valve with light membrane)

7717 Rx Provox® FreeHands HME Replacement Device
(Valve with medium membrane)

7721 Rx Provox® FreeHands HME Replacement Device
(Valve with strong membrane)

7730 Rx Provox® HME Cap (1 pc)

7251 Regular Round Adhesive (20 pcs)

7252 Regular Oval Adhesive (20 pcs)

7253 FlexiDerm Round Adhesive (20 pcs)

7254 FlexiDerm Oval Adhesive (20 pcs)

7255 OptiDerm Round Adhesive (20 pcs)

7256 OptiDerm Oval Adhesive (20 pcs)

7265 XtraBase Adhesive, Round (20 pcs)

7624 Rx LaryTube 8/36 with Ring (incl. 5 pcs HiFlow + 1 pc
FlexiDerm Oval)

7625 Rx LaryTube 8/55 with Ring (incl. 5 pcs HiFlow + 1 pc
FlexiDerm Oval)

7626 Rx LaryTube 9/36 with Ring (incl. 5 pcs HiFlow + 1 pc
FlexiDerm Oval)

7627 Rx LaryTube 9/55 with Ring (incl. 5 pcs HiFlow + 1 pc
FlexiDerm Oval)

7628 Rx	LaryTube 10/36 with Ring (incl. 5 pcs HiFlow + 1 pc FlexiDerm Oval)
7629 Rx	LaryTube 10/55 with Ring (incl. 5 pcs HiFlow + 1 pc FlexiDerm Oval)
7630 Rx	LaryTube 12/36 with Ring (incl. 5 pcs HiFlow + 1 pc FlexiDerm Oval)
7631 Rx	LaryTube 12/55 with Ring (incl. 5 pcs HiFlow + 1 pc FlexiDerm Oval)
7244	Provox® Cleaning Towel (200 pcs)

59420425	SkinPrep™ (50 pcs)
720400	SkinPrep™ (50 pcs) U.S. Version
59403125	Remove™ (50 pcs)
403100	Remove™ (50 pcs) U.S. Version











Provox® HME and Provox® FreeHands HME® are protected by US patents (5.042.468 and 5.738.095) and other patents and patents pending. Provox® XtraBase has a patent.










Provox®, FreeHands HME® and XtraBase® are registered trademarks owned by Atos Medical AB, Sweden.

Skin-Prep and Remove are trademarks of Smith & Nephew United, Inc., USA

CAUTION:

United States federal law restricts devices identified herein with the symbol "Rx" to sale, distribution and use by or on the order of a physician or licensed health care provider. The availability of these devices without prescription outside the United States may vary from country to country.

Speech valve membranes REF 7713, -14, -15

Disposable Provox FreeHands HME Cassettes REF 7712

Membrane Replacement Forceps REF 7719

Provox FreeHands speech valve unit REF 7710

Cleaning and Storage Box REF 7718

Safety Screwdriver REF 7719

Provox XtraBase Adhesive REF 7265

© Atos Medical AB, 200704A, REF 90750

ATOS

MEDICAL

Atos Medical AB
P.O. Box 183, SE-242 22 Hörby, Sweden
Tel: +46 (0)415 198 00 • Fax: +46 (0)415 198 98
info@atosmedical.com • www.atosmedical.com